

Powerful in the application –
soft to the tooth.

Main features
and advantages:

- Innovative ultrasonic technology for an optimal treatment result
- Continuously adjustable via KaVo foot control
- Highly precise, effective tooth cleaning controlled by linear oscillations
- 3 different Scaler tips
- Easy to use



KaVo. Dental Excellence.

Available tips: Prophylaxis



Piezo Scaler Tip 201

for the removal of supragingival calculus in all quadrants
1.007.4024



Piezo Scaler Tip 202

for the removal of supra- and subgingival calculus in all quadrants, especially in the approximal and sulcus areas
1.007.4026



Piezo Scaler Tip 203

for the removal of subgingival deposits on the root surfaces. Also suitable for periodontal preventative treatments
1.007.4028



Delivery contents:

- KaVo PIEZO*soft* hand piece with instrument tubing
- 1x PIEZO Scaler Tip 201
- Torque wrench
- Wear indicator card

Technical data

Drive type	piezoelectric
Oscillation form	linear
Oscillation frequency	24-32 kHz
Handpiece Ø	max. 17 mm
Handpiece length (without tip):	103 mm
Weight	50 g
Removal capacity	max. 8W



KaVo. Dental Excellence.